

DETRO MATIC ENZYME ENZYMATIC CLEANING SOLUTION FOR INSTRUMENT WASHER MACHINES, MANUAL, ULTRASONIC MACHINES

DESCRIPTION OF PRODUCT:

Enzymatic cleaning solution for automatic instrument washer machines, manual and ultrasonic machines

SPECIFICATION OF PRODUCTS:

Chemical-Physical Characteristic

Parameters	M.U	Specifications	Method of Analysis
Physical form	-	Liquid	DETC01
Appearance	-	Limpid	DETC03
Colour	-	Colourless	DETC05
Odour	-	Characteristic	DETC07
pH	pH unit, 20°C	8-11	DETC09
Density (20°C)	g/cm ³ , 20°C	1,05 ± 0,1 g/cm ³	DETC11
Solubility in water	-	Soluble	DETC13

*Confirmed by Detrox Quality Control Department.

DETRO MATIC ENZYME ENZYMATIC CLEANING SOLUTION FOR INSTRUMENT WASHER MACHINES, MANUAL, ULTRASONIC MACHINES

INTENDED USE OF PRODUCT:

Product Properties

- Compatible For All Washer Machines
- Compatible For All Water Hardness
- Cleaning and Disinfection of Heat-Resistant Instruments at 93 °C
- Cleaning of Thermolabile Instruments at 50-60 °C
- Protection Against Corrosion
- Low Usage Concentration
- Powerful Effect

Usage Areas

- It is used for machine enzymatic cleaning of medical instruments made from metal, glass, porcelain, plastic and other equipments.

Usage form

- At machine, it is used at % 0,2 - 0,8 (20 - 80 ml product to 10 lt water) concentration,
Manually it is used % 1 (10 ml product to 10 lt water) concentration.

Ingredients

- Enzymes (amylase, protease, lipase), amphoteric tensides, complexing agents, corrosion inhibitors

First aid informations:

- Eyes; In the case of eye contact, wash plenty of water around 5 min.
- Skin: Wash your skin immediately.
- Ingestion: Drink water and apply your doctor.

Netice TONBİL
General Manager

DETRO HEALTHCARE KİMYA SAN. A.Ş.
Atatürk Mahallesi Akşirler Menfesiler, Çarşıbaşı
No: 7 Etiler - Beşiktaş / İstanbul / Türkiye
Tel: 0212 659 77 62 - 0212 659 77 63
Faks: 0212 659 77 63 - 0212 659 77 63