

## DETRO ENZYM – INSTRUMENT&ENDOSCOPE CLEANER

### ***DESCRIPTION OF PRODUCT:***

Enzymatic cleaning solution for precleaning of instruments and endoscopes before disinfection process.

### ***SPECIFICATION OF PRODUCTS:***

#### Chemical-Physical Characteristic

Parameters	M.U	Specifications	Method of Analysis
Physical form	-	Liquid	DETC01
Appearance	-	Limpid	DETC03
Colour	-	Yellowish	DETC05
Odour	-	Characteristic	DETC07
pH	pH unit, 20°C	6,0-8,0	DETC09
Density (20°C)	g/cm <sup>3</sup> , 20°C	1,037 ± 0,1 g/cm <sup>3</sup>	DETC11
Solubility in water	-	Soluble	DETC13

\*Confirmed by Detrox Quality Control Department.

