

# DETRO MATIC NS NEUTRALIZATION SOLUTION FOR INSTRUMENT WASHER MACHINES

## ***DESCRIPTION OF PRODUCT:***

Neutralization solution for automatic instrument washer machine. It is used to neutralize instruments after alkaline cleaning in machines.

## ***SPECIFICATION OF PRODUCTS:***

### Chemical-Physical Characteristic

Parameters	M.U	Specifications	Method of Analysis
Physical form	-	Liquid	DETC01
Appearance	-	Limpid	DETC03
Colour	-	Colourless	DETC05
Odour	-	Characteristic	DETC07
pH	pH unit, 20°C	1,5-3	DETC09
Density (20°C)	g/cm <sup>3</sup> , 20°C	1,48 ± 0,1 g/cm <sup>3</sup>	DETC11
Solubility in water	-	Soluble	DETC13

\*Confirmed by Detrox Quality Control Department.

# DETRO MATIC NS NEUTRALIZATION SOLUTION FOR INSTRUMENT WASHER MACHINES

## **INTENDED USE OF PRODUCT:**

### Product Properties

- Removes Alkali Remnants
- Provides Effective Cleaning
- Compatible For All Washer Machines
- Compatible for All Water Hardness
- Protection Against Corrosion
- Low Usage Concentration
- Powerful Effect

### Usage Areas

- Detro Matic NS is used for removing alkali remnants after alkali cleaning of instruments at washer machines.

### Dosing:

- It is used at %0,05-0,1 concentration. (5-10 ml product into per 10 lt water)

### Ingredients

- Citric acid, complexing agents, corrosion inhibitors

### First aid informations:

- Eyes; In the case of eye contact, wash plenty of water around 5 min.
- Skin: Wash your skin immediately.
- Ingestion: Drink water and apply your doctor.

Netice TONBİL  
General Manager

DETRO HEALTHCARE KİMYA SAN. A.Ş.  
Atatürk Mahallesi / Ankara / Marmara Caddesi  
No: 7 Etiler / Beşiktaş / İstanbul / T.C.  
Tic. Sic. No: 271115 / Mersis No: 3401010000000007770001  
Fax: 0 212 659 77 62 - 0 212 659 77 63

**DETRO HEALTHCARE KİMYA SANAYİ A.Ş.**

☎ 212 659 77 62 📠 212 659 77 63

✉ info@detrox.com.tr - www.detrox.com.tr