

DETRO MATIC NF NEUTRALIZATION SOLUTION FOR INSTRUMENT WASHER MACHINES

DESCRIPTION OF PRODUCT:

Neutralization solution for automatic instrument washer machine. It is used to neutralize instruments after alkaline cleaning in machines.

SPECIFICATION OF PRODUCTS:

Chemical-Physical Characteristic

Parameters	M.U	Specifications
Physical form	-	Liquid
Appearance	-	Limpid
Colour	-	Colourless
Odour	-	Characteristic
pH	pH unit, 20°C	1,5-3
Density (20°C)	g/cm ³ , 20°C	1,48 ± 0,1 g/cm ³
Solubility in water	-	Soluble

*Confirmed by Detrox Quality Control Department.

DETRO MATIC NF NEUTRALIZATION SOLUTION FOR INSTRUMENT WASHER MACHINES

INTENDED USE OF PRODUCT:

Product Properties

- Removes Alkali Remnants
- Provides Effective Cleaning
- Compatible For All Washer Machines
- Compatible for All Water Hardness
- Protection Against Corrosion
- Low Usage Concentration
- Powerful Effect

Usage Areas

- Detro Matic NF is used for removing alkali remnants after alkali cleaning of instruments at washer machines.

Dosing:

- It is used at %0,05-0,1 concentration. (5-10 ml product into per 10 lt water)

Ingredients

- Phosphoric acid, complexing agents, corrosion inhibitors

First aid informations:

- Eyes; In the case of eye contact, wash plenty of water around 5 min.
- Skin: Wash your skin immediately.
- Ingestion: Drink water and apply your doctor.

Netice TONBİL
General Manager

DETROX HEALTHCARE KİMYA SANAYİ A.Ş.
Ankara Merkezli Tıbbi Aletler Üretimi
ve Tedavisi İçin Çözümler Üreten
Tıbbi Aletler Üretimi ve Tedavisi İçin Çözümler Üreten